Umbilical cord clamp

Clamp para cordão umbilical

- FR** Clamp de Bar ombilical
- DE** Nabelschnurklemme
- NL** Navelstrengklem
- IT** Morsetto per cordone ombelicale
- ES** Pinza para cordón umbilical
- SV** Navelsträngsklämma
- BG** Клампа за пъпна връв
- HR** Stezaljka za pupčanu vrpцу
- CS** Svorka na pupěčník
- DA** Navlestrengsklemmer
- ET** Nabanõõri klamber
- FI** Napanuoran puristin
- EL** Λαβίδα ομφάλιου λώρου
- HU** Köldökzsinórcsipesz
- LV** Nabassaites spailes
- LT** Virkštelės spaustukas
- PL** Zacisk do pępowiny
- RO** Clemă cordon ombilical
- SK** Svorka na pupočnú šnúru
- SL** Sponka za popkrovino

REF 489-051 STERILE EO CE 1639




Bastos Viegas, s.a.


Avenida da Fábrica, 298. 4560-164 Guilhufe, Penafiel - Portugal
Tel.: 255 729 500 - geral@bastosviegas.com - www.bastosviegas.com

IFU551999105, Rev.1 (25-10-2023)

1 <https://www.bastosviegas.com/downloads/ifu/IFUE551999105.pdf>

2

EN Umbilical cord clamp

INTENDED PURPOSE, INDICATIONS AND CLINICAL BENEFIT

Umbilical cord clamp is intended to clamp the umbilical cord after delivery to prevent bleeding from the cord stump after the cord had been cut.

INSTRUCTIONS FOR USE

- Check the thickness of the umbilical cord. In case it is a thin cord, it may be necessary to fold it before placing the clamp to obtain a bigger resistance to the clamp closing.
- Place the cord within the clamp, in the clamping zone (figure 1a).
- Check that there is not any object blocking the locking site (e.g. swabs, surgical instruments, etc.).
- Close the clamp until it locks and clicks (figure 2). Make sure that the clamp is well closed and that it is impossible to open it with the hands. Cut the umbilical cord.
- Confirm that there is no bleeding at the end of the cord. Check regularly for signs of bleeding due to a misplaced clamp or due to a damaged cord.

WARNINGS AND PRECAUTIONS

- Umbilical cord clamp is to be used exclusively by healthcare professionals.
- Always check the integrity of the locking system before use. Do not use if incomplete or blocked. Incorrect placement of the clamp on the umbilical cord may lead to an inadequate seal and consequent blood loss that might prove fatal for some babies.
- Do not re-use. Intended to be used one single time to prevent risks from contamination, including possible infections, and/or functional damage.
- Visually inspect the sterile barrier system prior to use. Do not use the device if you detect any breach in packaging integrity.
- Visually inspect the device. Do not use the device if you detect any sign of deterioration or defect that may compromise the performance.
- Dispose according to local laws rules, regulations and/or infection prevention standards for urban waste, to minimize potential risk of contamination, infection or pollution.

The occurrence of serious incidents is not expected with the use of the device if used correctly, but according to Regulation (EU) 2017/745 we must inform that serious incidents related with the use of the device, that are life-threatening or compromise health, including illness or permanent impairment of a body structure, must be notified to Bastos Viegas and to the Competent Authority of your Country.

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