EN Sterilization wraps

4501 - Sterilization wrap made of nonwoven SMS polypropylene, blue

4502 - Sterilization wrap made of nonwoven SMS polypropylene, green (alcohol repellent)

4503 - Sterilization wraps made of nonwoven SMS polypropylene, blue and green interleaved

450 - Creped medical grade paper (white, blue or green)

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450 - Papel crepado de grau médico (branco, azul ou verde)

EN STERILIZATION WRAPS INTENDED PURPOSE AND CLINICAL BENEFIT

Sterilization wraps are used for packaging of medical devices for sterilization in order to maintain the sterility of terminally sterilized medical devices until the point of use, in order to prevent patient infection.

Nonwoven SMS polypropylene sterilization wraps are designed for sterilization by Moist Heat (suitable for 121°C or 134°C; process in conformance with EN ISO 17665), ethylene oxide (process in conformance with EN ISO 11135) and vaporized hydrogen peroxide - VH2O2 (process in conformance with ISO 22441).

Crepe paper is designed for sterilization by Moist Heat (suitable for 121°C or 134°C; process in conformance with EN ISO 17665), ethylene oxide (process in conformance with EN ISO 11135) and radiation (process in conformance with EN ISO 11137).

PRODUCT INFORMATION

In compliance with the requirements from EN ISO 11607-1/-2 and EN 868-2, aiming to ensure microbial barrier properties and the subsequent sterility maintenance.

Designed for an easy and versatile aseptic presentation of the sterilized medical devices.

INSTRUCTIONS FOR USE

1-Visually inspect the sterilization wraps before use, in order to confirm there are no defects or foreign elements.

2-Select the size of the sterilization wrap according to size, shape and weight of the medical devices to be wrapped or based on guidelines within the organization.

The size should be selected to achieve adequate coverage of the articles being wrapped but it should not be so big that it has to be wrapped several times around the medical devices, as this could difficult the penetration of the sterilant agent.

3-Perform the proper wrapping technique, in order to allow aseptic presentation of the sterile content. The most common examples of wrapping techniques are the envelope fold and the square fold.

When two layers of sterilization wraps are used, the simultaneous double-wrapping technique (that refers to wrapping two layers of wrap material together simultaneously), or the sequential wrapping technique (that refers to wrapping two layers of wrap material individually) can be performed.

Note that it is essential to wrap the medical devices securely to prevent gaps, billowing and air pockets, but do not wrap too tightly as this could create holes or tears in the wrap. Furthermore, the sterilization wrap should be large enough to accommodate movement of the wrap during the sterilization cycle without ripping or tearing.

4- Close the wrapped packages with an appropriate closure system. The most common closure system is the indicator tape.

Note that the closure system should maintain the integrity of the package.

Do not use closures systems that avoid the penetration of the sterilant agent nor that compresses the package (e.g., ropes, strings, elastic bands, paperclips, staples or similar items).

5-Submit the package to the appropriate validated sterilization method.

6-After sterilization, visually inspect the package to ensure it is intact and the product is clean and dry. Wet packages are considered non-sterile.

7-Handle and storage the product carefully in order to maintain the integrity of the packaging.

If sterilization is carried out by a subcontracted unit, it is recommended that the packages are protected from possible contamination with a sterility maintenance cover.

INSTRUCTIONS FOR STERILIZATION

The below Moist Heat and EO sterilization process parameters have been validated according to EN ISO 17665, for Moist heat and EN ISO 11135, for Ethylene Oxide (EO).

MOIST HEAT:

CYCLE PARAMETERS	SPECIFICATION	DURATION
Vacuum/Air removal	≤100mbar x 3 cycles	-
Exposure Temperature	135-137°C	5 min
	OR	
	121-124°C	20 min

ETHYLENE OXIDE:

CYCLE PARAMETERS	SPECIFICATION
Pressure	45 - 1025 mbar
Humidity	95%HR
Injection rate	170mbar/min
Vacuum rate	50mbar/min
Exposure Temperature	55°C
Exposure time	1335 min.
Exposure EO concentration	750mg/L

The instructions provided above have been validated by Bastos Viegas. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result.

WARNINGS AND PRECAUTIONS

- Intended to be used exclusively by sterile processing personnel.
- Articles to be packed must be completely clean and dry.
- Avoid sharp edges and other surfaces that might disrupt the package. If it is not possible to avoid them, they should be protected.
- Keep the sterilization wraps in a clean and dry environment.
- We decline any responsibility in case of:
 - Innapropriate use;
 - Innapropriate storage or handling conditions;
 - Innapropriate sterilization conditions or sterilization method.
- Do not re-use. Intended to be used one single time to prevent risks from contamination, including possible infections, and/or functional damage.
- Do not resterilize.
- Visually inspect the device. Do not use the device if you detect any sign of deterioration or defect that may compromise the performance.
- Sterilization wraps maintain the product sterility up to 180 days.
- Dispose according to applicable legislation and/or infection prevention standards for healthcare waste, to minimize potential risk of contamination, infection or pollution.

The occurrence of serious incidents is not expected with the use of the device if used correctly, but according to Regulation (EU) 2017/745 we must inform that serious incidents related with the use of the device, that are life-threatening or compromise health, including illness or permanent impairment of a body structure, must be notified to Bastos Viegas and to the Competent Authority of your Country.

LABELLING SYMBOLS

1)UKCA marking 2)United Kingdom authorized representative

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