Bastos Vegas

Mêche de gazeTamponadebindeGaaswiekStriscia di garzaTira de gasaTamponad av gasvávМарля на лентаKrastaka gazaGázová páska

EN RIBBON GAUZE

With X-ray detectable thread / Without X-ray detectable thread
INTENDED PURPOSE AND CLINICAL BENEFIT
Ribbon gauze is intended to be used:
For packing natural body orifices and surgical wounds.
Mode of action includes absorption of drainage and fluids.
The application of pressure to achieve haemostasis or as a carrier for administration of medications.

INDICATIONS

-For nasal packing (treatment of epistaxis; for administration of medicinal substances).
-For ear packing (for administration of medicinal substances e.g. treatment of otitis externa).
-For packing of natural body orifices (e.g. for absorption of body fluids; for post-surgical packing; for stabilization of bony pieces).

INSTRUCTIONS FOR STERILIZATION

The device must be packaged within a sterile barrier system (SBS) permeable and compatible with Moist Heat and/or Ethylene Oxide (EO) gas sterilization and with the sterilization parameters that are specified in the following table. The system shall comply with EN ISO 16067-1/2 and CEN ISO 15135:2021. It is the responsibility of the processor to select an appropriate SBS and carry out the packaging under controlled environmental conditions. The below Moist Heat and EO sterilization process parameters have been validated according to EN ISO 16067-1/2 and EN ISO 15135, for EO.

Parameters for heat sterilization:
-Maximum density of the sterilization load must be 10 kg/m³.
-Packaging comprising a maximum of 3 levels, being paper/film peel-pack the SBS (cart. to EN ISO 16067-1/2 and EN 868-5), card and carton box the protective packaging.
-Sterility assurance level of 10⁻⁶.

CYCLE PARAMETERS	SPECIFICATION	DURATION
Vacuum /Air removal	≤100mbar x 3 cycles	-
Exposure temperature	135-137°C	5min.
	OR	
	121-124°C	20min.

Parameters for Ethylene Oxide sterilization:
-Maximum density of the sterilization load must be 150kg/m³ (in euro pallets).
-Packaging comprising a maximum of 3 levels, being paper/film peel-pack the SBS (cart. to EN ISO 16067-1/2 and EN 868-5), card and carton box the protective packaging.
-Sterility assurance level of 10⁻⁶.

CYCLE PARAMETERS	SPECIFICATION	DURATION
Minimum temperature during the preconditioning	20 °C	-
Temperature	33 - 46 °C	-
Humidity	50 ± 15 %RH	22-27h

CYCLE PARAMETERS	SPECIFICATION	DURATION
Vacuum /Air removal	≤100mbar x 3 cycles	-
Humidification/Dwell	60 - 95 %RH	60min.
Gas injection	350 mbar of N ₂ 480 mbar of EO	60min.
Exposure pressure	930 mbar	-
Exposure EO concentration	510 - 750 mg/L	173-1335min.
Exposure temperature	48 - 55 °C	-
EO removal	70mbar x 3 cycles	17min.
Flushing	70mbar x 3 Cycles	-

CYCLE PARAMETERS	SPECIFICATION	DURATION
Minimum temperature during the preconditioning	20 °C	-
Temperature	33 - 46 °C	-
Humidity	50 ± 15 %RH	22-27h

Two cycles, for both methods with the above parameters have been validated. The second sterilization cycle is only allowed for use after the first cycle. After the second cycle, processing, the device's shelf life will be limited to 5-years. This 5-year duration corresponds to the shelf life between the sterilized component and the validity of the non-sterile component.

Moist heat and EO sterilization processing shall be performed by trained personnel and by utilizing adequate and validated equipment's. The processor must comply with EN ISO 17665 (for Moist heat sterilization) and EN ISO 11135 (for EO sterilization) and any national legislation where the processing must occur. If transport is applicable, it must be validated acc. to international applicable standards (e.g., ASTM4169, ISTA).

INSTRUCTIONS FOR USE

After sterilizing:
-Assess the appropriateness of ribbon gauze for the medical intervention.
-If suitable, follow healthcare guidelines or your organization's protocols for being (ongoing monitoring, and changing) ribbon gauze.
-Use sterile scissors to cut the appropriate length of ribbon gauze.
-Maintain proper (including aseptic) technique to ensure favorable clinical outcome and patient safety, while preventing infection.
-If cutting the ribbon gauze, any open and unused ribbon gauze shall be discarded.
-Document the procedure in the patient's medical record, including the type and amount of ribbon gauze used and any observations or concerns of the ongoing care.

WARNINGS AND PRECAUTIONS
-Intended to be used sterile.
-Intended to be used exclusively by healthcare professionals.
-Visually inspect the device. If you use the device and you detect any sign of deterioration or defect that may compromise the performance of the device, do not use it.
-Do not re-use. Intended to be used one single time to prevent risks from contamination, including possible infections, and/or functional damage.
-Dispose according to applicable legislation and/or infection prevention standards for healthcare waste, to minimize potential risk of contamination, infection or pollution.

The occurrence of serious incidents is not expected with the use of the device if used correctly, but according to Regulation (EU) 2017/745 we inform that serious incidents may occur, such as:
-A life-threatening or compromise health, including illness or permanent impairment of a body part, if the device is not used according to the instructions for use.
-A life-threatening or compromise health, including illness or permanent impairment of a body part, if the device is not used according to the instructions for use.

LABELLING SYMBOLS
1) United Kingdom authorized representative 2) Nr of threads/cacm 3) Number of units 4) With X-ray detectable thread

Bastos Vegas,s.a.
Avenida da Fábrica, 298,
4560-164 Guilhufe, Penafiel- Portugal
Tel.: 255 729 100 - geral@bastosvegas.com
www.bastosvegas.com

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