
RETRACTOR ROUX


- FR Écarteur de Roux
- DE Roux-Haken
- NL Roux-retractor
- IT Divaricatore Roux
- ES Separador de Roux
- SV Roux-sårhake
- BG Петрактор Roux
- HR Retraktor Roux
- CS Rouxův retractor

RETRACTOR WEITLANER

- FR Écarteur de Weitlaner
- DE Weitlaner-Haken
- NL Weitlaner-retractor
- IT Divaricatore Weitlaner
- ES Separador de Weitlaner
- SV Weitlaner-sårhake
- BG Петрактор Weitlaner
- HR Retraktor Weitlaner
- CS Weitlanerův retractor

https://www.bastosviegas.com/downloads/ifu/UFU51999122.pdf

Bastos Viegas,s.a.
Avenida da Fábrica, 298,
4560-164 Guilhufe, Penafiel- Portugal
Tel: 255 729 500 - geral@bastosviegas.com
www.bastosviegas.com

UFU51999122. Rev. 1 (03-12-2025)

ROUX RETRACTOR WEITLANER RETRACTOR

INTENDED PURPOSE AND INDICATIONS

To hold back soft tissues and organs in order to expose the surgical field.

INSTRUCTIONS FOR STERILIZATION

Most heat (MH) and ethylene oxide (EO) sterilization must be conducted by trained personnel using validated protocols and sterilizers. The processor shall comply with ISO 17665-1/2/3 (for MH sterilization) and ISO 11335 (for EO sterilization), as well as any applicable national regulations. If transportation is involved, it must be validated according to internationally recognized standards (eg, ISTA, ASTM1469).

Devices must be packaged within a sterile barrier system (SBS) suitable for the chosen sterilization method and compatible with the respective process parameters. The packaging system must be designed to ensure the integrity of the device during and after sterilization, ensuring sterility to be achieved during the sterilization process and maintained throughout the product's shelf life.

The processor is responsible for selecting and validating an appropriate SBS for the specific sterilization method used. The packaging must be carried out under controlled environmental conditions to maintain product integrity.

During handling and use, the processor shall maintain environmental control to prevent contamination.

The sterilization process parameters for MH and EO, outlined below, have been validated in accordance with ISO 17665-1/2/3 (for MH) and ISO 11335 (for EO). However, the customer is responsible for validating their own sterilization process to ensure compliance with regulatory and operational requirements, as well as achieving a SAL of 10⁻⁶, considering their chosen SBS and the packaging levels.

Validated parameters for moist heat (MH) sterilization:

- The maximum allowable sterilization load density is 110 kg/m³.
- Maximum bioburden level of the article is 100 CFU/g.
- Packaging must consist of a minimum of two levels.
- First level: A validated paper-plastic SBS.
- Second level: A protective layer designed to preserve the bioburden level of the first level (SBS), such as a cardboard box or other suitable outer packaging.

Validated parameters for ethylene oxide (EO) sterilization:

- The maximum allowable sterilization load density is 180 kg/m³.
- Maximum bioburden level of the article is 100 CFU/g.
- Packaging must consist of a minimum of two levels and an outer protective layer designed to preserve the bioburden level of the first level (SBS).
- Third level: An outer protective packaging, such as a carton box, providing additional protection during handling and transport.

CYCLE PARAMETERS	SPECIFICATION	DURATION
Vacuum /Air removal	≤100mbar x 3 cycles	5*10min**
Exposure temperature	121-134°C	OR 20*40min**
Bubbling	≤100mbar	20*40min**

*If the minimum exposure time is used, a second sterilization cycle is permitted. Please refer to the information below for details on double sterilization.

Validated parameters for ethylene oxide (EO) sterilization:

- The maximum allowable sterilization load density is 180 kg/m³.
- Maximum bioburden level of the article is 100 CFU/g.
- Packaging must consist of a minimum of two levels and an outer protective layer designed to preserve the bioburden level of the first level (SBS).
- Third level: An outer protective packaging, such as a carton box, providing additional protection during handling and transport.

CYCLE PARAMETERS	SPECIFICATION	DURATION
Minimum temperature during the Preconditioning	33 – 46 °C	-
Relative Humidity	50 ± 5 %RH	22 - 72h
Vacuum /Air removal	≤100mbar x 3 cycles	60 min.
Humidification/Dwell	60 - 95 %RH	52 min.
Gas injection	480 mbar of EO	52 min.
Exposure pressure	930 mbar	-
Exposure EO concentration	580 - 750 mg/L	420-1335min.
Exposure temperature	48 - 55 °C	-
EO removal	100 mbar x 3 cycles	17 min.
Flushing	70mbar x 3 cycles	-

Temperature 28 - 42 °C Min. 11h

Speed of air extraction 10 - 20 m/s

The validated sterilization process permits up to two cycles using the specified parameters, with a restriction on moist heat sterilization as outlined in the *. A second sterilization cycle is only allowed for unused devices.

Post-sterilization, the device's shelf life is determined by the duration between the end of the product's intrinsic shelf life and the shelf life of the selected SBS. Additionally, the processor is responsible for ensuring compliance with the residual levels defined for the article, as specified in ISO 10993-7. Regarding storage conditions, all Bastos Viegas articles must be protected from sunlight and kept away from rain throughout their shelf life.

WARNINGS AND PRECAUTIONS

- Intended to be used sterile.
- Not intended to be used on hard tissues.
- Visually inspect the device. Do not use the device if you detect any sign of deterioration or defect that may compromise the performance.
- Intended to be used exclusively by healthcare professionals. Do not re-use. Intended to be used one single time to prevent risks from contamination, including possible infections, and/or damage.
- Dispose according to applicable legislation and/or infection prevention standards for healthcare waste, to minimize potential risk of contamination, infection or pollution.

The occurrence of serious incidents is not expected with the use of the device if used correctly, but according to Regulation (EU) 2017/475 as amended, the manufacturer must monitor the use of the device if there are life-threatening or compromise health, including lethal or permanent damage, or if there are serious risks to the health of patients.

The Competent Authority of your Country.

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