Ribbon gauze

Mecha de gaze

FR Mèche de gaze
DE Tamponadebinde
NL Gaaswiek
IT Striscia di garza
ES Tira de gasa
SV Tamponad av gasväv
BG Марля на лента
HR Trakasta gaza
CS Gázová páska

<https://www.bastosviegas.com/downloads/ifu/IFUe551999112.pdf>

Bastos Viegas,s.a.

Avenida da Fábrica, 298. 4560-164 Guilhufe, Penafiel- Portugal
Tel.: 255 729 500 - geral@bastosviegas.com
www.bastosviegas.com

IFU551999112. Rev. 3 (04-02-2025)

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EN RIBBON GAUZE

-WITH X-RAY DETECTABLE THREAD / -WITHOUT X-RAY DETECTABLE THREAD

INTENDED PURPOSE AND CLINICAL BENEFIT

Ribbon gauze is intended to be used:
For packing natural body orifices and surgical wounds.
Mode of action includes absorption of drainage and fluids, support, application of pressure to achieve haemostasis or as a carrier for administration of medicines.

INDICATIONS

-For nasal packing (treatment of epistaxis; for administration of medicinal substances).
-For ear packing (for administration of medicinal substances e.g. treatment of otitis externa).
-For packing of surgical wounds (e.g. for absorption of body fluids; for post-surgical packing; for stabilization of bony pieces).

INSTRUCTIONS FOR USE

-Assess the appropriateness of ribbon gauze for the medical intervention.
-If suitable, follow healthcare guidelines or your organization's protocols for using (including monitoring, and changing) ribbon gauze.
-If necessary, use sterile scissors to cut the appropriate length of ribbon gauze.
-Maintain proper (including aseptic) technique to ensure favorable clinical outcome and patient safety, while preventing infection.
-If cutting the ribbon gauze, any open and unused ribbon gauze shall be discarded.
-Document the procedure in the patient's medical record, including the type and amount of ribbon gauze used and any observations or instructions for ongoing care.

WARNINGS AND PRECAUTIONS

-Intended to be used exclusively by healthcare professionals.
-Visually inspect the sterile barrier system prior to use. Do not use the device if you detect any breach in packaging integrity.
-Visually inspect the device. Do not use the device if you detect any sign of deterioration or defect that may compromise the performance.
-Do not re-use. Intended to be used one single time to prevent risks from contamination, including possible infections, and/or functional damage.
-Dispose according to applicable legislation and/or infection prevention standards for healthcare waste, to minimize potential risk of contamination, infection or pollution.

The occurrence of serious incidents is not expected with the use of the device if used correctly, but according to Regulation (EU) 2017/745 we must inform that serious incidents related with the use of the device, that are life-threatening or compromise health, including illness or permanent impairment of a body structure, must be notified to Bastos Viegas and to the Competent Authority of your Country.

LABELLING SYMBOLS

1) Paper (identification of packaging material) 2) Opening direction 3) Nr of threads/sq cm 4) Number of units 5) Trademark 6)With X-ray detectable thread 7) United Kingdom authorized representative

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