# Skin staple remover

## Removedor de agrafos

FR Ôte-agrafes cutanées

DE Hautklammerentferner

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IT Rimuovi graffe

ES Quitagrapas

SV Agraffborttagare

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
## EN SKIN STAPLE REMOVER

### INTENDED PURPOSE AND INDICATIONS

Skin staple remover is intended for removal of surgical skin staples.

### INSTRUCTIONS FOR USE

- Inspect the staples and make sure they are perpendicular to skin before removal.
- Place the lower jaw of staple remover between staple and skin (figure 1).
- Squeeze and lift the staple remover perpendicularly to the skin.
- The staple will be bent open by squeezing (figure 2).
- Lift the open staple straight up from the tissue and discard according to local requirements.

### WARNINGS AND PRECAUTIONS

- Intended to be used exclusively by healthcare professionals.
- Visually inspect the sterile barrier system prior to use. Do not use the device if you detect any breach in packaging integrity.
- Visually inspect the device. Do not use the device if you detect any sign of deterioration or defect that may compromise the performance.
- Do not re-use. Intended to be used one single time to prevent risks from contamination, including possible infections, and/or functional damage.
- Dispose according to applicable legislation and/or infection prevention standards for healthcare waste, to minimize potential risk of contamination, infection or pollution.

The occurrence of serious incidents is not expected with the use of the device if used correctly, but according to Regulation (EU) 2017/745 we must inform that serious incidents related with the use of the device, that are life-threatening or compromise health, including illness or permanent impairment of a body structure, must be notified to Bastos Viegas and to the Competent Authority of your Country.

### LABELLING SYMBOLS

1)Trademark 2)Opening direction 3)United Kingdom authorized representative

4)Polyethylene (identification of packaging material)

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