Roux Retractor

Afastador Roux

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DE Roux-Haken

NL Roux-retractor

IT Divaricatore Roux

ES Separador de Roux

SV Roux-särhake

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
EN ROUX RETRATOR

WEITLANER RETRATOR

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Moist heat (MH) and ethylene oxide (EO) sterilization must be conducted by trained personnel using validated processes and sterilizers. The processor shall comply with ISO 17665-1/2/3 (for MH sterilization) and/or ISO 11135 (for EO sterilization), as well as any applicable national regulations. If transportation is involved, it must be done according to internationally recognized standards (e.g., ISTA; ASTM4169).

Devices must be packaged within a sterile barrier system (SBS) suitable for the chosen sterilization method and compatible with the respective process parameters. The chosen system must be validated according to ISO 1607-1/2 and CEN/ISO TS 16775, allowing sterility to be maintained throughout the sterilization process and maintained throughout the product's shelf life.

The processor is responsible for selecting and validating appropriate sterilization methods and methods used. The packaging must be tested under controlled environmental conditions to maintain product integrity.

During handling of the article, the processor shall maintain environmental control to prevent contamination.

The sterilization process parameters for MH and EO, and the maximum allowable sterilization load density with ISO 17665-1/2/3 (for MH) and ISO 11135 (for EO). However, the customer is responsible for validating their own sterilization process to ensure compliance with regulatory and operational requirements, as well as achieving a SAL of 10⁻⁶, considering their chosen SBS and the packaging levels.

Validated parameters for moist heat (MH) sterilization:

• The maximum allowable sterilization load density is 110 kg/m³.

• Maximum bioburden level of the article is 100 CFU/g.

• Packaging must consist of a minimum of two levels and a first level: A validated paper-plastic SBS.

• Second level: A protective layer designed to preserve the bioburden level of the first level (SBS), such as a cardboard box or other suitable outer packaging.

CYCLE PARAMETERS **SPECIFICATION** **DURATION**

Vacuum / Air removal ≤100mbar x 3 cycles -

Exposure temperature 135-137°C 5*-10min*

Drying OR 121-124°C 20*-40min*

 ≤100mbar ≥5min

*If the minimum exposure time is used, a second sterilization cycle is permitted. Please refer to the information below for details on double sterilization.

Validated parameters for ethylene oxide (EO) sterilization:

• The maximum allowable sterilization load density is 180 kg/m³.

• Maximum bioburden level of the article is 100 CFU/g.

• Packaging must consist of a minimum of two levels and a first level: A validated paper-plastic SBS.

• Second level: A protective layer designed to preserve the bioburden level of the first level (SBS).

• Third level: An outer protective packaging, such as a carbon box, providing additional protection during handling and transportation.

CYCLE PARAMETERS **SPECIFICATION** **DURATION**

PRECONDITIONING

Minimum temperature 21,4 °C -

Temperature 33 - 46 °C 22 - 72h

Relative Humidity 50 ± 15 %RH

STERILIZATION

Vacuum / Air removal -

Humidification / Dwell 60 - 95 %RH 60 min.

Gas injection 350 mbar of N₂O 52 min.

Exposure pressure 350 mbar of EO

Exposure temperature 48 - 55 °C

EO removal 100 mbar 17 min.

Flushing 70mbar x 3cycles

AERATION

Temperature 28 - 42 °C

Speed of air extraction 10 - 20 m/s

 Min. 11h

The validated sterilization process permits up to two cycles using the specified parameters, with a restriction on moist heat sterilization as outlined in the*. A second sterilization cycle is only allowed if the unit used does not exceed the shelf life of the selected SBS. Additionally, the processor is responsible for ensuring compliance with the residual levels defined for the article, as specified in ISO 10937-1. In all cases, all Bastos Viegas articles must be protected from sunlight and kept away from rain throughout their shelf life.

*If the minimum exposure time is used, a second sterilization cycle is permitted. Please refer to the information above on double sterilization.

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