



# OR Dressings

Material de cobertura cirúrgica

## LABELLING INFORMATION

Some of our products with double packing have an double control tag with the necessary information for an easy and safe use in the specific procedure in which they are to be used. Each tag is pre-cutted in two parts. Once apart, one tag goes along with the product to the operating room. The other remains out of the operating room for reconciliation of number of pieces at end of the operation

## INFORMAÇÃO SOBRE ROTULAGEM

*Alguns dos nossos produtos com dupla embalagem, tem etiqueta dupla de controlo com a informação necessária para uma fácil e segura utilização nos procedimentos específicos nos quais vão ser utilizados. Cada etiqueta é picotada em duas partes. Uma vez destacadas, uma parte vai junto com o produto para o bloco operatório. Outra parte da etiqueta fica fora do bloco operatório para reconciliação do numero de unidades no final da operação.*

Check number icon  
 Ícone verifique o número

Product drawing  
 Desenho do produto

Size & other product informations  
 Tamanho & outras informações do produto

Brand  
 Marca

Number of pieces  
 Número de unidades

Identification TAG number  
 Número de identificação da etiqueta

Pre-cutted  
 Pré-cortado

Datamatrix  
 Datamatrix

Reconciliation  
 Reconciliação

10 x 60 cm  
 8 ☞ 17 gx5pcs  
 Nr H 589080

10 x 60 cm  
 8 ☞ 17 gx5pcs  
 NrH 589080

10 x 60 cm  
 8 ☞ 17 gx5pcs  
 Nr H 589080